

Mouse Lptn (Lptn) Antibody Pair (STJAP00970)

GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Goat/Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Concentration	0.5 mg/mL/0.5-4ug/mL
Conjugation	Unconjugated/Biotin
Dilution Range	1:1000-1:10000
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	16963
Gene Symbol	Xcl1
UniProt ID	XCL1_MOUSE

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---