

Anti-Araucaria virus Nucleoprotein antibody (STJHVA00008)

GENERAL INFORMATION

Product Type	Primary Antibodies
Applications	ELISA/IHC/WB
Host / Source	Rabbit
Reactivity	Andes orthohantavirus

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Purified by antigen affinity column.
Dilution Range	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Isotype	IgG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--