

Anti-BCCV Gn/Glycoprotein N antibody (STJHVA00011)

GENERAL INFORMATION

Product Type	Primary Antibodies
Applications	ELISA/IHC/WB
Host / Source	Rabbit
Reactivity	Black Creek Canal orthohantavirus/BCCV/Black Creek Canal virus

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Purified by antigen affinity column.
Dilution Range	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Isotype	IgG
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--