

Anti-Hantaan virus/HTNV GP antibody [3G1] (STJHVA00012)

GENERAL INFORMATION

Product Type	Primary Antibodies
Applications	ELISA/FCM/IF/WB
Host / Source	Human
Reactivity	Hantaan virus/strain 76-118/Korean hemorrhagic fever virus

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Dilution Range	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--