

Anti-Hantaan virus/HTNV Glycoprotein N/Gn antibody [AH100] (STJHVA00014)

GENERAL INFORMATION

Product Type	Primary Antibodies
Applications	ELISA/IF/NEUT/IP
Host / Source	Human
Reactivity	Hantaan virus

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Dilution Range	ELISA: 1:2000-1:20000
Formulation	0.01M PBS pH 7.4
Isotype	IgG
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--