

## Anti-SEOV Gc antibody (STJHVA00026)

### GENERAL INFORMATION

Product Type	Primary Antibodies
Applications	ELISA/IHC/WB
Host / Source	Rabbit
Reactivity	Seoul virus/strain R22

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Purified by antigen affinity column.
Dilution Range	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Isotype	IgG
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--