

## Anti-botA/BOTOX nanobody [SAA1450] (STJN000001)

### GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA/IHC/WB
Host / Source	Alpaca
Reactivity	C.botulinum

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.05 mg/mL
Conjugation	Unconjugated
Purification	Purified by Nickel column.
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-His
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

