

## Anti-Mal f 1/Allergen Pit o 1 nanobody [SAA0639] (STJN000005)

### GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA
Host / Source	Alpaca
Reactivity	M.furfur

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.98 mg/mL
Conjugation	Unconjugated
Purification	Purified by Nickel column.
Isotype	VHH-8His-Cys-tag
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

