

## Anti-CNR1/CB1 nanobody [SAA1274] (STJN000170)

### GENERAL INFORMATION

|               |               |
|---------------|---------------|
| Product Type  | Nanobodies    |
| Applications  | ELISA/RIA/SPR |
| Host / Source | Alpaca        |
| Reactivity    | Human         |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 1 mg/mL  |
| Conjugation         | Unconjugated   |
| Purification        | Purified by Nickel column.   |
| Formulation         | 0.01M PBS, pH 7.4.   |
| Isotype             | VHH-8His-Cys-tag   |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">1268</a>       |
| Gene Symbol | <a href="#">CNR1</a>       |
| UniProt ID  | <a href="#">CNR1_HUMAN</a> |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

