

## Anti-CD184/CXCR4 nanobody [SAA1262] (STJN000240)

### GENERAL INFORMATION

|               |            |
|---------------|------------|
| Product Type  | Nanobodies |
| Applications  | ELISA/RIA  |
| Host / Source | Alpaca     |
| Reactivity    | Human      |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                   |
| Concentration       | 1 mg/mL                                                                                                                                                      |
| Conjugation         | Unconjugated                                                                                                                                                 |
| Purification        | Purified by Nickel column.                                                                                                                                   |
| Formulation         | 0.01M PBS, pH 7.4.                                                                                                                                           |
| Isotype             | VHH-8His-Cys-tag                                                                                                                                             |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">7852</a>        |
| Gene Symbol | <a href="#">CXCR4</a>       |
| UniProt ID  | <a href="#">CXCR4_HUMAN</a> |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

