

## Anti-CD114/CSF3R nanobody [SAA1247] (STJN000295)

### GENERAL INFORMATION

|               |              |
|---------------|--------------|
| Product Type  | Nanobodies   |
| Applications  | ELISA/IF/SPR |
| Host / Source | Alpaca       |
| Reactivity    | Human        |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 1 mg/mL  |
| Conjugation         | Unconjugated   |
| Purification        | Purified by Nickel column.   |
| Formulation         | 0.01M PBS, pH 7.4.   |
| Isotype             | VHH-8His-Cys-tag   |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">1441</a>        |
| Gene Symbol | <a href="#">CSF3R</a>       |
| UniProt ID  | <a href="#">CSF3R_HUMAN</a> |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

