

Anti-ANGPTL4 nanobody [SAA1258] (STJN000300)

GENERAL INFORMATION

| | |
|---------------|------------|
| Product Type | Nanobodies |
| Applications | ELISA |
| Host / Source | Alpaca |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | Purified by Nickel column. |
| Formulation | 0.01M PBS, pH 7.4. |
| Isotype | VHH-8His-Cys-tag |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 51129 |
| Gene Symbol | ANGPTL4 |
| UniProt ID | ANGL4_HUMAN |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

