

Anti-HRSV F/Fusion glycoprotein F0-strain A2 nanobody [SAA0982] (STJN000423)

GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA
Host / Source	Alpaca
Reactivity	HRSV-A

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by Nickel column.
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-8His-Cys-tag
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

