

Anti-Clostridium botulinum botA/BOTOX-strain Hall/ATCC 3502/NCTC 13319/Type A nanobody [SAA0922] (STJN000460)

GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA
Host / Source	Alpaca
Reactivity	C.botulinum

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by Nickel column.
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-8His-Cys-tag
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

