

Anti-Vibrio cholerae serotype O1 higB-2-strain ATCC 39315/EI Tor Inaba N16961 nanobody [SAA0877] (STJN000486)

GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA
Host / Source	Alpaca
Reactivity	V.cholerae serotype O1

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by Nickel column.
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-8His-Cys-tag
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

