

Anti-Influenza A virus NP/Nucleoprotein protein-strain A/Wilson-Smith/1933 H1N1 nanobody [SAA1425] (STJN000532)

GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA/FC/IF
Host / Source	Alpaca
Reactivity	Influenza A

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.93 mg/mL
Conjugation	Unconjugated
Purification	Purified by Nickel column.
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-8His-Cys-tag
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

