

Anti-Human CD332/FGFR2 (IIIb) nanobody [SAA2047] (STJN000562)

GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA/FC
Host / Source	Alpaca
Reactivity	Human/Cynomolgus

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by Nickel column.
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-8His-Cys-tag
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

TARGET INFORMATION

Gene ID	2263
Gene Symbol	FGFR2
UniProt ID	FGFR2_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

