

## Anti-Human CD16a/FCGR3A nanobody [SAA2051] (STJN000566)

### GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA/FC
Host / Source	Alpaca
Reactivity	Human/Cynomolgus

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by Nickel column.
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-8His-Cys-tag
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

### TARGET INFORMATION

Gene ID	<a href="#">2214</a>
Gene Symbol	<a href="#">FCGR3A</a>
UniProt ID	<a href="#">FCG3A_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

