

Anti-CD206/MRC1 nanobody [3.49#] (STJN000570)

GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA
Host / Source	Alpaca
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by Nickel column.
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-8His-Cys-tag
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

TARGET INFORMATION

Gene ID	4360
Gene Symbol	MRC1
UniProt ID	MRC1_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

