

Anti-VHH nanobody [cAb-Lys2] (STJN000638)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/Isotype Control
Host / Source	Alpaca

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Isotype	VHH-8His-Cys-tag
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at-20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--