

Anti-HSV gD/US6 nanobody [SAA2253] (STJN000652)

GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA/SPR/IF
Host / Source	Alpaca
Reactivity	HSV

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Purified by Nickel column.
Dilution Range	ELISA: 1:2000-1:20000
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--