

## Anti-MERS-CoV-RBD nanobody [VHH-55] (STJN000655)

### GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA/Neutralization
Host / Source	Alpaca
Reactivity	MERS-CoV

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Purified by Nickel column.
Dilution Range	ELISA: 1:2000-1:20000
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-His
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--