

Anti-IL23R nanobody [SAA2269] (STJN000664)

GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA/FCM/SPR
Host / Source	Alpaca
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Dilution Range	ELISA: 1:2000-1:20000
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-Fc
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	149233
Gene Symbol	IL23R
UniProt ID	IL23R_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--