

## Anti-Nipah virus/NiV Fusion glycoprotein nanobody [DS90] (STJN000671)

### GENERAL INFORMATION

|               |             |
|---------------|-------------|
| Product Type  | Nanobodies  |
| Applications  | ELISA       |
| Host / Source | Alpaca      |
| Reactivity    | Nipah virus |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                           |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                |
| Concentration       | Lot specific                                                                                                                                                              |
| Conjugation         | Unconjugated                                                                                                                                                              |
| Purification        | Purified by Nickel column.                                                                                                                                                |
| Dilution Range      | ELISA: 1:2000-1:20000                                                                                                                                                     |
| Formulation         | 0.01M PBS, pH 7.4.                                                                                                                                                        |
| Isotype             | VHH-His                                                                                                                                                                   |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |

### ADDITIONAL INFORMATION

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|