

## Anti-CD158a/KIR2DL1 nanobody [SAA2281] (STJN000674)

### GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA
Host / Source	Alpaca
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Purified by Nickel column.
Dilution Range	ELISA: 1:2000-1:20000
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-His
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">3802</a>
Gene Symbol	<a href="#">KIR2DL1</a>
UniProt ID	<a href="#">K12L1_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--