

Anti-Mu-type opioid receptor/OPRM1/Mu OR nanobody [SAA2285] (STJN000676)**GENERAL INFORMATION**

Product Type	Nanobodies
Applications	ELISA/IF/FCM
Host / Source	Human
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Dilution Range	ELISA: 1:2000-1:20000
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-hFc
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	18390
Gene Symbol	Oprm1
UniProt ID	OPRM_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--