

Anti-Gal-1/LGALS1 nanobody [SAA2310] (STJN000687)**GENERAL INFORMATION**

Product Type	Nanobodies
Applications	ELISA/WB
Host / Source	Alpaca
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Purified by Nickel column.
Dilution Range	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Formulation	0.01M PBS pH 7.4
Isotype	VHH-His
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3956
Gene Symbol	LGALS1
UniProt ID	LEG1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--