

## Anti-FZD3/Frizzled-3 nanobody [SAA2319] (STJN000689)

### GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA
Host / Source	Human
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Dilution Range	ELISA: 1:2000-1:20000
Formulation	0.01M PBS pH 7.4
Isotype	VHH-Fc
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">7976</a>
Gene Symbol	<a href="#">FZD3</a>
UniProt ID	<a href="#">FZD3_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--