

## Anti-VHH nanobody [cAb-Lys2] (STJN000691)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/Isotype Control
Host / Source	Alpaca

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.1 mg/mL
Conjugation	Unconjugated
Purification	Purified by Nickel column
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--------------------------------------------------------------------------------------------------------------------