

Anti-VEGF/ANG-2 nanobody [BI 836880] (STJN000692)

GENERAL INFORMATION

Product Type	Nanobodies
Applications	ELISA/WB
Host / Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Purified by ion exchange chromatography.
Dilution Range	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-VHH-VHH
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--