

## Anti-SV2A nanobody [SAA2262] (STJN000706)

### GENERAL INFORMATION

Product Type	Nanobodies
Applications	IP
Host / Source	Alpaca
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Purified by Nickel column.
Formulation	0.01M PBS, pH 7.4.
Isotype	Alpaca
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">64051</a>
Gene Symbol	<a href="#">Sv2a</a>
UniProt ID	<a href="#">SV2A_MOUSE</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--