

ACE2 protein (C-His&Avi) (STJP000026)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4. or Supplied as a 0.22 Mu m filtered solution in PBS, pH 7.4.
Molecular Weight	86.2kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	59272
Gene Symbol	ACE2
UniProt ID	ACE2_HUMAN
Immunogen Region	Gln18-Ser740
Immunogen Sequence	QSTIEEQAKTFLDKFNHEAE DLFYQSSLASWNYNTNITEE NVQNMNAGDKWSAFLKEQS TLAQMYPLQEIQLNLTVKLQL QALQQNGSSVLSSEDKSKRLN TILNTMSTIYSTGKVCNPDN PQECLLLEPGLNEIMANSLD YNERLWAWESWRSEVVGKQLR PLYEEYVVLKNEMARANHYE DYGDYWRGDYEVNGVDGYDY SRGQLIEDVEHTFEEIKPLY EHLHAYVRAKLMNAYPSYI

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

