

MANF protein (STJP000139)

GENERAL INFORMATION

Product Type	Proteins
--------------	----------

PRODUCT PROPERTIES

Concentration	1 mg/ml
---------------	---------

Formulation	PBS pH 7.4, concentration
-------------	---------------------------

Purification	Purified by ion-exchange chromatography and gel-filtration from serum-free CHO growth media
--------------	---

Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles
---------------------	--

TARGET INFORMATION

Gene ID	7873
---------	----------------------

Gene Symbol	MANF
-------------	----------------------

UniProt ID	MANF_HUMAN
------------	----------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---