

2B4 protein (STJP000141)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	CHO cells

PRODUCT PROPERTIES

Formulation	PBS pH 6.0
Purification	Purified by Ni-affinity chromatography, ion-exchange chromatography and desalting chromatography from serum-free CHO growth media
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID	18106
Gene Symbol	Cd244
UniProt ID	CD244_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---