

CD48 protein (STJP000142)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	CHO cells

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	Purified by ion-exchange chromatography, hydrophobic interaction chromatography and gel-filtration from serum-free CHO growth media
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID	12506
Gene Symbol	Cd48
UniProt ID	CD48_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---