

GFRA-1 protein (STJP000149)

GENERAL INFORMATION

Product Type Proteins

PRODUCT PROPERTIES

Concentration 0.9 mg/ml

Formulation PBS pH 7.4

Purification Purified by Ni-affinity chromatography and gel-filtration from serum-free CHO growth media

Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID [2674](#)

Gene Symbol [GFRA1](#)

UniProt ID [GFRA1_HUMAN](#)

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.