

GFRA-2 protein (STJP000153)

GENERAL INFORMATION

Product Type	Proteins
--------------	----------

PRODUCT PROPERTIES

Concentration	1.3 mg/ml
---------------	-----------

Formulation	PBS pH 7.4
-------------	------------

Purification	Purified by Ni-affinity chromatography and gel-filtration from serum-free CHO growth media
--------------	--

Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles
---------------------	--

TARGET INFORMATION

Gene ID	2676
---------	----------------------

Gene Symbol	GFRA3
-------------	-----------------------

UniProt ID	GFRA3_HUMAN
------------	-----------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---