

## Cardiac Troponin complex protein (STJP000154)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	CHO cells

### PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	25 mM Tris, pH 7.5, 500 mM NaCl, 5 mM CaCl <sub>2</sub> , 15 mM Beta-mercaptoethanol
Purification	Antibody-affinity purification following desalting
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

### TARGET INFORMATION

Gene ID	<a href="#">7137</a>
Gene Symbol	<a href="#">TNNI3</a>
UniProt ID	<a href="#">TNNI3_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---