

ACE2 protein (hFc) (STJP000155)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Human

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	1. PBS pH 7.4, 2. 20 mM citrate, 150 mM NaCl; pH 6.0
Purification	Purified by Protein A affinity chromatography following gel filtration. Protein is sterile-filtrated through 0.22 μ m filter. NA
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID	59272
Gene Symbol	ACE2
UniProt ID	ACE2_HUMAN
Immunogen Region	18-740 aa

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---