

NGF protein (STJP000157)

GENERAL INFORMATION

Product Type Proteins

PRODUCT PROPERTIES

Concentration 1 mg/ml

Formulation PBS pH 7.4

Purification Purified by ion-exchange chromatography and gel-filtration from serum-free CHO growth media. Protein is sterile-filtrated through 0.22 μ m filter

Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID [4803](#)

Gene Symbol [NGF](#)

UniProt ID [NGF_HUMAN](#)

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.