

## Artemin protein (STJP000161)

### GENERAL INFORMATION

Product Type Proteins

### PRODUCT PROPERTIES

Concentration 1 mg/ml

Formulation 10 mM citrate; 150 mM NaCl pH 4.5

Purification Purified from serum-free CHO growth media by ion-exchange chromatography and gel-filtration. Protein is sterile-filtered through 0.22  $\mu$ m filter

Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

### TARGET INFORMATION

Gene ID [9048](#)

Gene Symbol [ARTN](#)

UniProt ID [ARTN\\_HUMAN](#)

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.