

## Recombinant SARS-CoV-2 Spike S1 protein (STJP000164)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	SARS-CoV-2

### PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	NA Metal-affinity chromatography following gel filtration. Protein is sterile-filtrated through 0.22 µm filter. NA
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

### TARGET INFORMATION

Immunogen Region	14-681 aa
------------------	-----------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---