

Recombinant SARS-CoV-2 Spike S2 ECD-His protein (STJP000165)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	SARS-CoV-2

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	NA Metal-affinity chromatography following gel filtration. Protein is sterile-filtrated through 0.22 μ m filter. NA
Storage Instruction	Store at -70°C upon receipt. It is strongly recommended to avoid repeated freeze-thaw cycles. Aliquot the protein, freeze quickly in liquid nitrogen and thaw only once for your experiments.

TARGET INFORMATION

Immunogen Region	693-1218 aa
------------------	-------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---