

SARS-CoV-2 Trimeric Spike protein (STJP000167)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	CHO cells

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	NA Purified by Ni-affinity chromatography and gel-filtration from serum-free CHO growth media, sterile filtrated NA
Storage Instruction	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---