

SARS-CoV-2 Spike S1 Indian B.1.617.2 protein (Delta) (STJP000171)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	CHO cells

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	Metal-affinity chromatography following gel filtration. Protein is sterile-filtrated through 0.22 μ m filter.
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---