

## Recombinant Interleukin 6 protein protein (STJP000183)

### GENERAL INFORMATION

Product Type Proteins

Host / Source CHO cells

### PRODUCT PROPERTIES

Concentration 1 mg/ml

Formulation PBS pH 7.4

Storage Instruction Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.