

## SARS-CoV-2 Trimeric Spike Californian B.1.429.1 protein (STJP000185)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	CHO cells

### PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	Purified by Ni-affinity chromatography and gel-filtration from serum-free CHO growth media, sterile filtrated
Storage Instruction	Store at-20 °C to-70 °C long term. Avoid multiple freeze-thaw cycles.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---