

Human CD32a (FcγRIIa) protein (Recombinant) {Biotin} (STJP000197)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised 0.2 µg/ml filtered PBS solution, pH7.2.
Molecular Weight	Recombinant CD32a is a monomer protein consisting of 192 amino acid residue subunits, due to glycosylation migrates as an approximately 25 kDa protein on SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	2212
Gene Symbol	FCGR2A
UniProt ID	FCG2A_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--