

Human EpCAM protein (Recombinant) {Biotin} (STJP000203)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant human EpCAM is a protein consisting of 249 amino acid residue subunits, due to glycosylation migrates as an approximately 35kDa protein on reduced SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	4072
Gene Symbol	EPCAM
UniProt ID	EPCAM_HUMAN
Immunogen Region	GLN24-LYS265, ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--