

Human Interleukin-15 protein (Recombinant) {Biotin} (STJP000207)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS solution pH7.5.
Molecular Weight	Native human Interleukin-15 is generated by the proteolytic removal of the signal peptide and propeptide with a theoretical molecular mass of approximately 13kDa. Recombinant Interleukin-15 Fc chimera is a homodimer protein and migrates an approximat
Storage Instruction	Can be stored in working aliquots at 4°C for one month, or at -20°C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	3600
Gene Symbol	IL15
UniProt ID	IL15_HUMAN
Immunogen Region	Mature chain

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--